

Parallel Chamber Extension

Design & Function

Our parallel chamber extension is designed to fit into the proximal luer connection of ModuFlo sets and contains everything required for safe secondary infusions from a spiked source:

- ⇒ Integral Infusion Chamber
- ⇒ Y-Connector
- ⇒ One-way check valve on chamber tubing
- ⇒ Precision roller valve and C-Clamp


The extension is fully compatible with infusion pumps, and offers a variety of placement options. Extension chamber height can be adjusted to suit your flow rate requirements.

It is your ideal way to connect two bags together and retain full control over infusion rates.

Uses

Uses include:

- ⇒ IV Medication by infusion
- ⇒ Secondary fluid infusion (eg plasma expander with saline)
- ⇒ Second bag for overnight volume requirements
- ⇒ Ideal connection to make up your own infusion set-up with dedicated IV access - eg for chemotherapy



Product Summary and Codes

Extension	Drop Size (drops/ml)	Length (cm)	Priming		Pliancy	Product Code
			ModuFlo side	FULL set		
Parallel Connector	20	75	1.7	4.0	High	ICYE20CV