

Injection Caps

- Male luer-lock cap with non-latex injection site
- Allows needled access into capped extensions or catheters
- Low profile and low priming volume
- Great value, especially in our bulk pack
- Ideal for capping catheters, extension lines & T-connectors and other open female luers that require injection access
- Available in packs of 10 or boxes of 100



Closed Caps

- Blind-ended caps with both male & female luer connections
- Connects to either female or male luers (or both)
- Excellent to cap dispensed liquids in syringes to clients
- Maintains sterility of luer connections, syringe hubs and catheter hubs
- Available in packs of 10 or boxes of 100



Swivel Connector

- Male & Female luer connector with silicone 'O' ring seal
- Allows 360° rotation either side of connector to prevent lines twisting
- Ideal for mobile animals that tend to twist giving sets
- Excellent for procedures with needled syringe that still require movement

Please ensure luer connections are firmly tightened before use



Double male luer connector

- Connector with male luers at both ends
- Connects two female luers together
- Packs of 5 individual sterile connectors



Cap for needle-free valve

- NON-activating cap for our needle-free valves
- Protects valve surface from gross contamination
- Ideal when exercising or bathing patient
- Available sterile or non-sterile



Product Summary and Codes

Cap or connector	Pack Size	Sterile?	Product Code
Injection Cap	10	Y	INJCAP10
Injection Cap - multi pack	100	Y	INJCAP-M
Closed Cap	10	Y	CLOSED10
Closed Cap - multi pack	100	Y	CLOSED-M
Swivel Connector	1	Y	CSW-MF-1
Double male luer connector	5	Y	C-MMLL-5
Cap for needle-free valve: Single, sterile	1	Y	FNVCAP-1
Cap for needle-free valve: Pack of 20, non-sterile	20	NO	FNVCAP20