

Mounting systems and pole stands

Mantis Mount (Registered Design: GB3020495)

- Lightweight, quality stainless steel clamp
- Mounts infusion pumps onto kennel doors
- New design features cross-brace for improved strength
- Leaves floor space free and uncluttered
- Riser rod available - fits into clamp to hold fluid bag
- Adaptor bar available for doors with vertical rails
- Fits standard (Shoreline, Geeling, Technik etc) & most custom kennels
- Custom versions to suit your pump and kennel are available. Made to order in powder coated mild steel
- When ordering with one of our pumps, just ask for the 'P' version to receive a significant saving over standard pricing.




Transfusion Stand mobile

- Stainless steel, double drip hook
- Height adjustable from 52" - 82" (132 - 208 cm)
- Base manufactured in mild steel, white powder coat finish
- Four free running castors



Infusion Stand mobile

- Economically priced
- Available with either two or four stainless steel hooks
- Base manufactured in solid aluminium, with 5 wheels for extra stability
- Height adjustable from 40" - 78" (100 - 195 cm)



Wall mounted Drip Hook

- White powder coated finish
- Supplied with all fittings



Product Summary and Codes

Holder	Hooks	Wheels	Product Code
Mantis Mount kennel door clamp	0	n/a	KC2004 (P)
Mantis Mount with riser rod holder for fluid bag	2	n/a	KC3004 (P)
Transfusion stand mobile	2	4	J2016
Infusion Stand mobile; TWO hooks	2	5	J2017
Infusion Stand mobile; FOUR hooks	4	5	J2017/4
Wall mounted Drip Hook	1	0	J2016W