

Jugular Catheters: Kits available

Basic Kits (-JC)

- Soft tipped polyurethane Seldinger catheter
- Guide wire (J & Straight tipped)
- Seldinger needle
- Bridge clamp (for securing part placed catheter with sutures)
- Scalpel blade and handle
- Dilator (for skin & vessel dilation)
- Non-Latex caps to close female luers



Product Summary & Codes

Custom-made catheters are available to suit your needs

Type	Lumens	O/d Gauge	Length (cm)	Channel gauges (needle-gauge)	Product Code
Single	1	3 Fr	8	20	1308-JC
Single	1	3 Fr	10	20	1310-JC
Single	1	3 Fr	13	20	1313-JC
Single	1	3 Fr	20	20	1320-JC
Single	1	4 Fr	10	18	1410-JC
Single	1	4 Fr	13	18	1413-JC
Single	1	5 Fr	15	16	1515-JC
Single	1	5 Fr	30	16	1530-JC
Single	1	6 Fr	20	14	1620-JC
Double	2	4 Fr	10	22/22	2410-JC
Double	2	4 Fr	12	22/22	2412-JC
Double	2	5 Fr	15	18/20	2515-JC
Double	2	7 Fr	20	16/16	2720-JC
Double	2	7 Fr	30	16/16	2730-JC
Triple	3	5.5 Fr	15	20/22/22	3515-JC
Triple	3	7 Fr	20	16/18/18	3720-JC
Triple	3	7 Fr	30	16/18/18	3730-JC

Needle-Free Kits (-NF)

- Contents as for basic kit as above, but with:
 - ⇒ **Sterile OpSite** dressing to protect catheter entry site
 - ⇒ **Needle-free valves** instead of injection caps, giving the best safety, sterility and ease of use



Product Summary & Codes

Type	Lumens	O/d Gauge	Length (cm)	Channel Gauges (needle-gauge)	Product Code
Single	1	3 Fr	10	20	1310-NF
Single	1	4 Fr	10	18	1410-NF
Single	1	5 Fr	15	17	1515-NF
Single	1	7 Fr	20	15	1720-NF
Double	2	4 Fr	12	16/20	2412-NF
Double	2	5 Fr	15	16/18	2515-NF
Double	2	7 Fr	20	14/14	2720-NF
Triple	3	5.5 Fr	15	18/22/22	3515-NF
Triple	3	7 Fr	20	16/18/18	3720-NF

CODES: Please note that our coding system has been simplified. All catheters have a 4 digit code followed by a two letter suffix. The first digit indicates the number of lumens; 2nd is the o/d Fr gauge and 3rd & 4th digits are the length (cm). The 'JC' suffix is a plain cap, whilst the 'NF' indicates a needle-free kit.