

Sundry Products


Delta Fix Catheter Dressing

- Sterile dressing for catheter sites
- With protective pad for improved comfort
- Can apply topical antiseptics if required
- Supplied individually wrapped in peel cover
- Available as singles or in boxes of 50



Lachrymal/Irrigating Cannula

- Stainless steel blunt-tipped irrigating cannula
- Angled profile for ease of use
- 27G - 14 G sizes to suit most needs
- Single sterile blister packed item



27 - 17 G are 60mm long
14 G is 70 mm long

Corrugated Drainage Sheet

- Corrugated PVC sheet suitable for non-cavity drainage
- PVC softens in contact with body warmth for patient comfort



Corrugated Drain Sheet with radiopaque marker

- Can be cut and stacked to regulate capillary action
- Available as single double pouched items

C-Clamps

- C-Clamps for chest drain closure
- Single sterile clamps in packs of 5
- Cost effective chest drain safety
 - ⇒ Medium: drains up to 18 Fr
 - ⇒ Large: drains 16 - 28 Fr



Cyanoacrylate Glue

- Ideal bonding agent for PVC drains
- Apply to all PVC joins for 100% security
- Handy 3g size for repeated applications
- Store in PVC bag in FRIDGE for maximum shelf life.



Pressure Cuff Infuser

- Cuff to pressurise I/V bags for rapid fluid administration
- Simple manometry bulb and release valve
- Built-in hook to hang from I/V pole
- Combine this with the **HiFlo** infusion line (page 7) for the most rapid fluid flow available in small animals - guaranteed!

- Pressure indicator & warning
- For 1,000 ml bags



Product Summary & Codes

Item	Colour	Size	Pack Size	Product Code
Delta Fix sterile dressing	n/a	n/a	1	DELTAFIX
Delta Fix sterile dressing - multi pack	n/a	n/a	50	DELTA-50
Lachrymal cannula: 27 G	purple	27G	1	LC4003.27
Lachrymal cannula: 22 G	black	22 G	1	LC4003.22
Lachrymal cannula: 19 G	brown	19 G	1	LC4003.19
Lachrymal cannula: 17 G	white	17 G	1	LC4003.17
Lachrymal cannula: 14 G	green	14 G	1	LC4003.14
Corrugated Drainage Sheet	clear/green	220 x 25 mm	1	DRS3059S
Super Glue	n/a	3 gm	1	J1706
C-Clamp; medium (up to 18 Fr)	red	medium	5	C-CLAMP-M5
C-Clamp; large (16 - 28 Fr)	white	large	5	C-CLAMP-L5
Pressure cuff infuser bag	white	1,000 ml	1	CUFF1000