


Sundry Products

Needle-free Injection Spike

- Use for sterile access into bag or bottle
- Needle-free valve is swabable, so does not need capping
- Spike is vented so can be used on non-collapsible containers
- Vent is filtered, so no risk of contaminating container
- Supplied as a single sterile unit in a blister pack



Mini dispensing pin

- Ideal for small bottles
- Filtered vent means no contamination of bottle content
- More compact than NFVSPIKE
- Fitted with cap, but ideally suits a needle-free valve for a totally sealed system content



Sterile isopropyl swabs

- Individual swabs containing isopropyl alcohol
- Foil wrapped and gamma sterilised
- Ideal for swabable needle-free valves and other IV access
- Perfect for immuno-compromised patients
- Supplied in boxes of 100 swabs



Syringes for syringe drivers

- 20 & 50/60ml luer-lock syringes for use with syringe drivers
- Smooth, 3-part, non-latex syringes
- Perfect accuracy with Pilot & Agilia range of syringe drivers
- Luer lock prevents accidental luer disconnection



PF5 syringes (20 & 60 ml)

Product Summary & Codes

Item	Capacity (ml)	Pack Size	Product Code
Needle-free injection spike	-	1	NFVSPIKE
Mini dispensing pin	-	3	MINIPIN3
Steret Alcohol swabs	-	100	STERETS
Spritze Syringe: 50 ml luer lock	50/60	1	SPRITZE
PF5 syringe, 20 ml, SINGLE	20	1	S-PF5-20
PF5 syringe, 20 ml, Pack of 12	20	12	SPF52012
PF5 syringe, 60 ml, SINGLE	50/60	1	S-PF5-60
PF5 syringe, 60 ml, Pack of 8	50/60	8	SPF56008