

Foley Catheters

Latex Foley Catheters (siliconised)

Features & Details

- Siliconised latex gives:
 - ⇒ Similar biocompatibility to silicone
 - ⇒ Economy of latex
- High performance valves for reliable two way action
- Soft and pliant for patient comfort
- Double-pouched for ease of use and sterility
- Small sizes (6-10Fr) have a PVC stiffening stylet to aid placement
- Available in a wide range of sizes
- High balloon volume available in 16 Fr catheter




Product Summary & Codes

Size	Stylet?	Length (cm)	Balloon Vol	Product Code
6 Fr	YES	25	1.5 ml	FC2006-06
8 Fr	YES	25	3-5 ml	FC2006-08
10 Fr	YES	25	3-5 ml	FC2006-10
12 Fr	No	34	5-15 ml	FC2006-12
16 Fr	No	34	30-50 ml	FC2006-16
20 Fr	No	34	5-15 ml	FC2006-20

100 % Silicone Foley Catheters

Features & Details

- 100 % silicone gives the best:
 - ⇒ biocompatibility
 - ⇒ placement duration
- Superbly economical silicone catheters
- 55 cm long; suitable for canine catheterisation
- Stiffening stylet supplied with all catheters
- Reliable two way valve action
- Supplied double pouched for simple, sterile placement
- Available in 6, 8 and 10 Fr catheter sizes; low balloon volumes
 Other sizes are available. We would welcome your enquiry



Product Summary & Codes

Size	Stylet?	Length (cm)		Balloon Vol	Product Code
		Total	Effective		
6 Fr	YES	55	52	1.5 ml	FC5506-S
8 Fr	YES	55	52	3 ml	FC5508-S
10 Fr	YES	55	52	3 ml	FC5510-S