

PLAIN Chest Drains: TROCHAR



Trochar drain (2012T)

Design

- DOUBLE pouched soft PVC drain
- Clear distance markings for accurate placement
- Radiopaque line for confirmation of placement on X-ray
- Sharp trochar for easy direct puncture
- Trochar-Drain interface ensures minimal tissue friction
- Standard connector allows attachment to your normal device,¹ or remove connector and use your preferred drain connector²
- Available as a single item
- An economical alternative when you need direct thoracic puncture

Details (2012T)



Trochar drain tip, drain holes and distal marker

- Wide double side holes for prolonged patency
- Connector already fitted to drain tubing
- Suitable for placement at routine thoractomy, or emergency thoracic puncture
- Contoured handle makes placement easy and comfortable
- Use as supplied, or use one of our needle-free connectors to give you full flexibility of connection, with peace of mind on security



Contoured handle for easy placement

Notes:

1. For secure attachment, we recommend you use cyanoacrylate glue on all connectors.
2. The needle-free drain connectors (DR-FNV-1, 2 & 3) allow a quick & easy connection to your preferred drain size. Please contact us for further details on the best drain connection for your needs.
3. 'Proximal' refers to length from proximal most fenestration hole

Product Summary and Codes

Drain Type	Gauge Fr	Length (cm)		Product Code
		Total	Proximal ³	
TROCHAR	12	16	13	DR2012T12
TROCHAR	14	16	13	DR2012T14
TROCHAR	16	18	15	DR2012T16
TROCHAR	18	22	18	DR2012T18
TROCHAR	20	32.5	28.5	DR2012T20
TROCHAR	24	32.5	28.5	DR2012T24
TROCHAR	28	32.5	28.5	DR2012T28