

Needle-Free KITS: TROCHAR



20 Fr Needle-free Trochar set (2012T20K)

Design

- DOUBLE pouched soft PVC drain
- Radiopaque line for confirmation of placement on X-ray
- Clear distance markings for accurate placement
- Wide double side holes for prolonged patency
- Available as a single item
- Drain material is the same as plain set (page 34)

DETAILS: TROCHAR set (2012T-K)

- Simple connection for secure assembly
- C-clamp fitted for security and attachment point
- Sharp trochar for easy thoracic puncture
- Suitable for placement by direct puncture
- Female luer connects to any luer-fitted device
- Valve connections depend on drain sizes:
 - ⇒ 12 & 14 Fr: Male & Female luer with connecting tube
 - ⇒ 16 Fr: Female luer
 - ⇒ 18+ Fr: Standard connector and DR-FNV-3
- Connection is needle-free, for easier aspiration and kernal management
 - ⇒ Valve closes automatically when detaching syringe or luer
 - ⇒ Self sealing means no need to cap the end to maintain sterile fluid path
 - ⇒ No risk of accidental opening



- Notes:
1. For secure attachment, it is recommended you use cyanoacrylate glue on trochar kit connection
 2. 'Proximal' refers to length from proximal most fenestration hole
 3. An Emergency Procedure Pack is available containing all items necessary for placing Trochar kits in an emergency - see page 37

Product Summary and Codes

Drain Type	Gauge Fr	Length (cm)		Pack Size	Product Code
		Total	Proximal ²		
TROCHAR	12	16	13	1	DR2012T12K
TROCHAR	14	16	13	1	DR2012T14K
TROCHAR	16	18	15	1	DR2012T16K
TROCHAR	18	22	21	1	DR2012T18K
TROCHAR	20	32.5	28.5	1	DR2012T20K
TROCHAR	24	32.5	28.5	1	DR2012T24K
TROCHAR	28	32.5	28.5	1	DR2012T28K