

Blood Collection, Transfer and Emergency Packs

Blood Collection Bags

- Your choice of 350 or 450 ml capacity
- Single, Twin or 10 pack for maximum convenience
- CPDA anti-coagulant for up to 28 day blood storage
- Bags are DOUBLE wrapped to ensure optimum anticoagulant protection
- FOIL outer wraps give maximum shelf life and do not require refrigeration
- Sealing clips available for clamping collection tube securely (packs of 20)
- Combine with our other transfusion products for improved transfusion efficiency (transfusion lines, needle-free spike for closed system bag access etc)



Presentations

- **450 ml bags:**
 - ⇒ **SINGLE** bag; outer wrap in **FOIL**
 - ⇒ **TWIN** bags; outer wrap in **FOIL**
- **350 ml bags:**
 - ⇒ **SINGLE** bag; outer wrap in **PVC**
 - ⇒ **TWIN** bags; outer wrap in **FOIL**
 - ⇒ **TEN-PACK**; outer wrap in **FOIL**

NOTE: • Prime collection tubing/syringes with anticoagulant prior to blood collection • BLOOD: ANTICOAGULANT volume ratio is **7.1 : 1** for syringe collections

Dry Transfer Bags

- Plain bags for use in blood separation
- 120 cm distance marked tubing
- Bag spike at end of tubing for blood bag access
- 150 and 300 ml bag sizes



Product	Capacity	Pack Size (Bags)	AntiCoag	Outer Wrap	Tube (cm)	Product Code
	Vol (ml)		Vol (ml)			
Blood Collection Bag	350 ml	1	49	PVC	150	T3031350
Blood Collection Bags	350 ml	2	49	FOIL	150	T3031-02
Blood Collection Bags	350 ml	10	49	FOIL	150	T3031-10
Blood Collection Bag	450 ml	1	63	FOIL	150	T3031.450.01
Blood Collection Bags	450 ml	2	63	FOIL	150	T3031.450.02
Dry Transfer Bag	150 ml	1	n/a	PVC	120	T3035.150
Dry Transfer Bag	300 ml	1	n/a	PVC	120	T3035.300
Clips for Blood Tube	n/a	20	n/a	PVC pack	n/a	T3030-20

Emergency Procedure Pack

- Contains all vital items you need for that urgent procedure:
 - ⇒ Outer tray drape (50 x 50 cm) wrapping all items together
 - ⇒ Inner drape (50 x 50 cm) with 6 cm diameter fenestration
 - ⇒ 3 Metric mono nylon, swaged 38mm TS
 - ⇒ Scalpel; # 11 blade; gauze swabs
 - ⇒ 20 ml syringe; 21 G needle
- Double-pouched pack for easy handling in an emergency
- Perfect for placing:
 - ⇒ Chest drain sets
 - ⇒ Tracheostomy tubes
 - ⇒ Jugular catheters, etc



Product Summary & Code

Product	Pack Size	Product Code
Emergency Procedure Pack	1	P-EPK-01