


OXYGEN DELIVERY & Feeding Tubes

Nasal Prong Set (twin bore)

- Twin bore sets for simple and economical oxygen supply
- Provides oxygen directly into nostrils, or use as 'flow by' supply
- Available with adult or paediatric prongs to suit most patient sizes
- Connects to cone connector (C-CONE-1)



Product Summary & Code

Product	Length (cm)	Size	Outlets	Pack Size	Product Code
Nasal Oxygen Prong Set	205	Adult	2	1	NP-2008A
Nasal Oxygen Prong Set	205	Paediatric	2	1	NP-2008P

Nasal Oxygen Catheters

- Single lumen catheters to supply oxygen directly to nasopharynx
- Two lateral holes and open end for good oxygen dispersal
- Moderate stiffness for easy placement
- Three sizes to suit most patient sizes
- Connect 6 & 8 Fr to Cone Connector (C-CONE-1); 10 Fr to a male luer




Product Summary & Code

Product	Length (cm)	Size (mm)	Outlets	Pack Size	Product Code
Nasal Oxygen CATHETER	37	6 Fr	3	1	NO201106
Nasal Oxygen CATHETER	39	8 Fr	3	1	NO201108
Nasal Oxygen CATHETER	39	10 Fr	3	1	NO201110

Nasal Tubes (feeding & oxygen)

- Soft, frosted PVC nasal tubes suitable for both feeding and oxygen delivery
- Stiffening stylet supplied for easy placement but maximum patient comfort
- Female luer lock for secure attachment and use
- Two lateral holes
- Other sizes and lengths available; please enquire
- Connects directly to male luer: - Oxygen Extension Set - NOE2100L
- Syringe



Product Summary & Code

Product	Length (cm)	Size (mm)	Outlets	Pack Size	Product Code
Nasal Feeding Tube	40	4 Fr	2	1	NG2038.04
Nasal Feeding Tube	50	6 Fr	2	1	NG2038.06
Nasal Feeding Tube	50	8 Fr	2	1	NG2038.08
Nasal Feeding Tube	100	10 Fr	2	1	NG2038.10
Nasal Feeding Tube	50	10 Fr	2	1	NG2038.10-5