

Silicone Feeding Tubes (OFTSXXXX)

Priming Volumes

Product Code	Description	Priming Volume (ml)
OFTS0835	Silicone feeding tube, 3 hole; 8 Fr, 35 cm	0.60
OFTS1040	Silicone feeding tube, 3 holes; 10 Fr, 40 cm	1.12
OFTS1225	Silicone oesophagostomy tube, 3 holes; 12 Fr, 25 cm	0.95
OFTS1250	Silicone feeding tube, 4 holes; 12 Fr, 50 cm	1.85
OFTS1425	Silicone oesophagostomy tube, 3 holes; 14 Fr, 25 cm	1.12
OFTS1460	Silicone feeding tube, 4 holes; 14 Fr, 60 cm	3.00
OFTS1630	Silicone oesophagostomy tube, 3 holes; 16 Fr, 30 cm	1.95
OFTS1660	Silicone feeding tube; 16 Fr, 60cm, 4 holes	4.05
OFTS2040	Silicone feeding tube, 4 holes; 20 Fr, 40 cm	4.20
OFTS2065	Silicone feeding tube, 4 holes; 20 Fr, 65 cm	7.60